Your Path Through the Sea

RBRsolo T

# Small Temperature Recorder

## Compact, accurate, dependable

The RBRsolo T is our most compact, lightweight, and versatile single channel recorder. The logger offers flexible measurement schedules, fast sampling up to 2Hz, large memory, extra power for extended deployments, fast USB download for large data files, and temperature accuracy of ±0.002°C.

#### **Features**

- Compact
- Lightweight
- High accuracy
- 30M measurements
- Up to 2Hz sampling
- Fast USB download speed



The RBRsolo T makes it easy to configure the optimum sampling regime for your measurements whether it is moored, towed, or profiling. The large data storage capacity and fast download ability facilitate long deployments with higher sampling rates. The RBRsolo T is small and streamlined making it the ideal instrument for shallow and deep water applications, down to 1700m (depth). The unique temperature sensor end cap facilitates flow at the sensor to reduce response time. Newly designed circuitry improves signal to noise, power consumption and allows for increased sampling rates. A dedicated desiccant holder makes it simple to use and replace desiccant before each deployment.

The RBRsolo T is one of the most flexible single channel recorders available from RBR and like all RBR instruments, the calibration coefficients are stored with the logger. Dataset export to Matlab™, Excel™, or text files makes post processing with your own algorithms simple and easy.







Your Path Through the Sea

# Small Temperature Recorder

## Compact, accurate, dependable

## Specifications

#### **Physical**

Power: 3.6V LiSOCI, AA cell

Communication: Fast USB

Storage: ~30M readings Clock Accuracy: ±60 seconds/year

Depth Rating: 1700m
Diameter: 25.4mm
Length: 240mm
Weight (air): 126g
Weight (water): 20g

#### **Temperature**

Range:  $-5^{\circ}\text{C}$  to  $35^{\circ}\text{C}$  Accuracy:  $\pm 0.002^{\circ}\text{C}$  Resolution:  $< 0.00005^{\circ}\text{C}$ 

Time Constant: 1.0s

Drift: 0.002°C/year

### **Deployment duration**

Sampling at 2Hz: 118 days (20.3M readings)
Sampling at 1s: 235 days (20.3M readings)
Sampling at 5s: 3 years (19.0M readings)

